

D396

Pad-Thermosol Range

Pad-Thermosol Range. Used for chemical impregnation or dyeing, squeezing and after treatment by infra-red preheater, intermediate dryer and thermosoling zone. Test fabrics 230mm wide by 400mm long are carried through the system by chains. Nip pressure, speed, heating capacity and processing temperatures are all variable to suit the fabrics being processed. Laboratory standard compressed air supply required.

WEIGHT	1600 kg	3520 lb		
DIMENSIONS	Width	Depth	Height	
	1200 mm 48 inch	3500 mm 140 inch	2900 mm 114 inch	
ORDERING Info.	D396	Pad-Thermosol Range		
	G265	Silent Laboratory Compressor (Optional) Larger capacities available on request		



SDL Atlas, Ltd

PO Box 162, Crown Royal, Shawcross St., Stockport SK1 3JW,
Tel: +44 (0) 161 480 8485 Fax: +44 (0) 161 480 8580
Internet: <http://www.sdlatlas.com>

SDL Atlas, LLC

3934 Airway Drive, Rock Hill, SC 29732, USA
Tel: +1 (803) 329-2110 Fax: +1 (803) 329-2133
Internet: <http://www.sdlatlas.com>

SDL Atlas, Ltd

1/F (South-East) & 2F, Shenjian Mansion, Central District
(West), Hi-Tech Park, Nanshan, Shenzhen, 518057, P.R.C.
Tel: +86 (755) 2671 1168 Fax: +86 (755) 2671 1337
Internet: <http://www.sdlatlas.com.cn>

SDL Atlas, Ltd

3J, Garment Centre, 576 Castle Peak Road, Kowloon, Hong
Kong
Tel: (852) 2173 8611 Fax: (852) 2370 2949
Internet: <http://www.sdlatlas.com>